



KenoseptTM L

SPRAY ON ALCOHOL BASIS FOR SURFACE TREATMENT

- ▶ easy and safe to use
- ▶ appropriate for all surfaces
- ▶ tested and proven effective against Influenza A/H1N1 (Mexican flu / swine flu)

Effective surface and hand disinfection

Kenosept L is the ideal solution in preventing the indirect transfer of micro-organisms via surfaces, contact points and hands.

Fast action

Kenosept L is -thanks to its unique composition- extremely appropriate for the fast treatment of tables, chairs and all washable surfaces.

Kenosept L has a bactericidal, fungicidal and virucidal action and complies to the norms EN 1040, EN 1276, EN 1500 and EN 1650.

Kenosept L contains alcohol and chlorhexidine which ensure a long lasting action.

Kenosept L is more than 90% biodegradable.

Kenosept L has been tested and proven effective against both the human as the swine strain of the Influenza virus A/H1N1. A copy of the report is available upon request.



Safe in use

Kenosept L is safe for all materials.
Kenosept L contains no fragrance.
Kenosept L is extremely soft for the hands.

Directions for use

Kenosept L is ready for use and can be applied in an (automatic) dispenser or with a dosing pump.

Treatment of surfaces:

- Spray the surfaces to treat with Kenosept L until completely covered
- Let the product dry
- Dry off with dry paper cloth if necessary

Hygienic hand disinfection:

- Apply 5 ml Kenosept L in the palm of the hand, rub hands well, also the space between the fingers
- Rub until completely dry

Chirurgical hand disinfection:

- Wash the hands with surgical soap in advance, rinse well with clear water
- Dry hands
- Apply 2 x 5 ml Kenosept L on hands and lower part of the arms, rub until completely dry

Packaging:

- 1 L
- 5 L

